




























Mode	Picture	Description	Price/USD
FS-610		High speed handpiece connector 2sets	
		Low speed handpiece connector 1set	
		3-way syringe (cold and hot) 1pcs	
		Operation lamp 1set	
		Film viewer 1set	
		Saliva ejector (strong and weak) 1set	
		Automatic cup filler 1set	
		rotatable ceramic spittoon 1set	
		Built-in water purified system 1set	
		24V D.C noiseless motor 1set	
		Sanitary cushion and backrest 1set	
		Doctor's stool 1pcs	
		FS-620A	
High speed handpiece connector			
Low speed handpiece connector			
3-way syringe (cold and hot)			
Operation lamp			
Film viewer with 3 memory			
Saliva ejector (strong and weak)			
Automatic cup filler			
rotatable glass spittoon			
Built-in water purified system			
24V D.C noiseless motor			
Sanitary cushion and backrest			
Doctor's stool			
FS-620		Motor dental chair control system with three memory	
		High speed handpiece connector	
		Low speed handpiece connector	
		3-way syringe (cold and hot)	
		Operation lamp	
		24V film viewer	
		Saliva ejector (strong and weak)	
		Automatic cup filler	
		Rotatable ceramic spittoon	
		Built-in water purified system	
		24V D.C noiseless motor	
		Sanitary cushion and backrest	
		Doctor's stool	
	\$2,007.00		

FS-620T		Motor dental chair control system with three memory		
		High speed handpiece connector		
		Low speed handpiece connector		
		3-way syringe (cold and hot)		
		Operation lamp		
		24V film viewer		1se
		Saliva ejector (strong and weak)		1set
		Automatic cup filler		1set
		Rotatable ceramic spittoon		1set
		Built-in water purified system		1set
		24V D.C noiseless motor		1set
		Sanitary cushion and backrest		
		Doctor's stool		
FS-630		Motor dental chair control system with three memory		
		High speed handpiece connector		
		Low speed handpiece connector		
		3-way syringe (cold and hot)		
		Operation lamp		
		24V film viewer		
		Saliva ejector (strong and weak)		
		Automatic cup filler		
		Rotatable glass spittoon		
		Built-in water purified system		
		Danish linak 24V noiseless DC motor		
		Sanitary cushion and backrest		
		Doctor's stool		1pcs
FS-630T		Motor dental chair control system with three memory		
		High speed handpiece connector		
		Low speed handpiece connector		
		3-way syringe (cold and hot)		
		Operation lamp		
		24V film viewer		
		Saliva ejector (strong and weak)		
		Automatic cup filler		
		Rotatable glass spittoon		
		Built-in water purified system		
		Danish linak 24V noiseless DC motor		
		Sanitary cushion and backrest		
		Doctor's stool		
FS-510		Arm rest		2set
		headrest		1pcs
		backrest		1set
		seatrest		1set
				\$775.00

		foot control	1set
		24V D.C noiseless motor	1set
FS-520		Arm rest	2set
		headrest	1pcs
		backrest	1set
		seatrest	1set
		foot control	1set
		24V D.C noiseless motor	1set
Opitional Item Price		up mounted type	196/set
		curinglight	69/set
		Woodepecker N2 detachable scaler	177/set
		oral camera syestem	584/set
		handpeice	196/set
		compressor	228/pcs
		fibre handpiece tube system	133/set
		Armrest	12.7/pcs
		leather cushion	260/set
Dental Oil-free Air Compressor			
FS-1.5EW-30		Voltage: 220V/110V 50Hz/60Hz Exhaust: 70L/min pressure: Noise: 56-65dB Supply for : 1 pcs dental unit	Frequency: Power : 550W Rated exhaust 0.8MPa Bottle volume : 30L
FS-1.5EW-30A (stainless steel tank)		Voltage: 220V/110V 50Hz/60Hz Exhaust: 70L/min pressure: Noise: 56-65dB Supply for : 1 pcs dental unit	Frequency: Power : 550W Rated exhaust 0.8MPa Bottle volume : 30L
FS-2EW-35		Voltage: 220V/110V Power : 840W Rated exhaust pressure: Noise: 56-65dB Supply for :	Frequency: 50Hz/60Hz Exhaust: 100L/min 0.8MPa Bottle volume : 35L 2 pcs dental unit
FS-2EW-35A (stainless steel tank)		Voltage: 220V/110V Power : 840W Rated exhaust pressure: Noise: 56-65dB Supply for :	Frequency: 50Hz/60Hz Exhaust: 100L/min 0.8MPa Bottle volume : 35L 2 pcs dental unit

FS-3EW-45		Voltage: 220V/110V Frequency: 50Hz/60Hz Power : 1100W Exhaust: 140L/min Rated exhaust pressure: 0.8MPa Noise: 56-65dB Bottle volume : 45L Supply for : 3 pcs dental unit	
FS-3EW-55		Voltage: 220V/110V Frequency: 50Hz/60Hz Power : 1100W Exhaust: 140L/min Rated exhaust pressure: 0.8MPa Noise: 56-65dB Bottle volume : 55L Supply for : 3 pcs dental unit	
FS-4EW-65		Voltage: 220V/110V Frequency: 50Hz/60Hz Power : 1680W Exhaust: 200L/min Rated exhaust pressure: 0.8MPa Noise: 56-65dB Bottle volume : 65L Supply for : 4 pcs dental unit	
FS-6EW-90		Voltage: 220V/110V Frequency: 50Hz/60Hz Power : 2520W Exhaust: 300L/min Rated exhaust pressure: 0.8MPa Noise: 56-65dB Bottle volume : 90L Supply for : 6 pcs dental unit	
FS-8EW-120		Voltage: 220V/110V Frequency: 50Hz/60Hz Power : 3316W Exhaust: 400L/min Rated exhaust pressure: 0.8MPa Noise: 56-65dB Bottle volume : 180L Supply for : 8 pcs dental unit	
FS-10EW-180		Voltage: 220V/110V Frequency: 50Hz/60Hz Power : 4200W Exhaust: 500L/min Rated exhaust pressure: 0.8MPa Noise: 56-65dB Bottle volume : 180L Supply for : 10 pcs dental unit	
FS-20EW-350		Voltage: 220V/110V Frequency: 50Hz/60Hz Power : 8400W Exhaust: 1000L/min Rated exhaust pressure: 0.8MPa Noise: 56-65dB Bottle volume : 350L Supply for : 20 pcs dental unit	
Dental Suction			
FS-550		Voltage: 220V/110V Frequency: 50Hz/60Hz Power : 550W vacuum pressure : 200mmHg Air flow: 1000L/Min Motor RPM:2800RPM	





FS-1100		Voltage: 220V/110V Frequency: 50Hz/60Hz Power : 840W vacuum pressure : 400mmHg Air flow: 2000L/Min Motor RPM:5600RPM	
FS-2EW-35/FS-550		Air compressor and suction integral Voltage: 220V/110V Frequency: 50Hz/60Hz Power : 550W vacuum pressure : 200mmHg Air flow: 1000L/Min Motor RPM:2800RPM	
Dental curing light			
FS-30		Halogen Curing light unit	
FS-060		LED wireless Curing light ,time setting :5s 10s 20s ,large capacity battery, a full charge can be used more than 430 times continuously underr 10s working time mode .power input:AC100V-240V 50/60HZ ,light output:850-1000mW/cm2	
FS-012A		LED wireless Curing light ,time setting :5s 10s 20s ,large capacity battery, a full charge can be used more than 430 times continuously underr 10s working time mode .power input:AC100V-240V 50/60HZ ,light output:850-1000mW/cm2	
FS-013		wired Curing light ,time setting :5s 10s 20s ,large capacity battery, a full charge can be used more than 430 times continuously underr 10s working time mode .power input:AC100V-240V 50/60HZ ,light output:850-1000mW/cm2	
FS-080		LED wireless Curing light ,time setting :5s 10s 20s ,large capacity battery, a full charge can be used more than 430 times continuously underr 10s working time mode .power input:AC100V-240V 50/60HZ ,light output:850-1000mW/cm2	

FS-090		LED wireless Curing light ,time setting :5s 10s 20s ,large capacity battery, a full charge can be used more than 430 times continuously underr 10s working time mode .power input:AC100V-240V 50/60HZ ,light output:850-1000mW/cm2
FS-022		Unifies both whitening accelerator and LED curing light Strong optical output power, maximum to 1800mW/cm2. Five treatment modes: Turbo, Normal, Soft start, Pulse, Pulse soft start Time setting: 0~60s for LED curing light, and 0~30 mins for whitening Full color LCD display Interchangeable battery pack, with both seat charger



Compressor

Mode	Picture	Name	FOB USD Price	Mode	Picture	Name	FOB USD Price
FS-MDT01		mobile dental unit	\$273.00	FS-PDT01		portable dental unit	
FS-550		oilless air compressor 550W motor	\$111.00	FS-840		oilless air compressor 840W motor	
FS-550		oilless air compressor 550W motor	\$123.00	FS-840		oilless air compressor 840W motor	

Steam Sterilizers

Mode	Picture	Name	FOB USD Price	Mode	Picture	Name	FOB USD Price
FS-8		Steam Sterilizers 8L	\$1,163.00	FS-12		Steam Sterilizers 12L	
FS-18		Steam Sterilizers 18L	\$1,239.00	FS-23		Steam Sterilizers 23L	
















Whitening unit

Mode	Picture	Name	FOB USD Price	Mode	Picture	Name	FOB USD Price
FS-WA05		whitening unit	\$470.00	FS-WA06		whitening unit	

FS-WA07		whitening unit	\$465.00		FS-WA08		whitening unit	\$374.00
X-ray								
Mode	Picture	Name	FOB USD Price		Mode	Picture	Name	FOB USD Price
FS-R0A		Portable X-ray	\$1,163.00		FS-R08-1		X-ray wall type	
FS-R0B		Portable X-ray (with battery)	\$1,356.00		FS-R08		X-ray	
FS-R05		automatic film processing	\$413.00		FS-RS01		digital dental x-ray	
FS-R03		Manual film processing	\$118.00		FS-RD01		x-ray film reader	
Dental air prophy								
Mode	Picture	Name	FOB USD Price		Mode	Picture	Name	FOB USD Price
FS-SG01		dental air prophy	\$38.80		FS-SG02		dental air prophy	
FS-F086		Lemon Sand Power(130g)	\$1.65		FS-F087		Mint Sand power (130g)	
Ultrasonic scaler								
Mode	Picture	Name	FOB USD Price		Mode	Picture	Name	FOB USD Price
FS-G3		ultrasonic scaler	\$103.00		UDS-K		ultrasonic scaler	
UDS-J		ultrasonic scaler	\$142.00		N2		ultrasonic scaler	
FS-AC01		Air scler	\$83.00					
Dental stool								
Mode	Picture	Name	FOB USD Price		Mode	Picture	Name	FOB USD Price

FS-F042		Small light 220V	\$36.20		FS-F366		Light 220V	\$77.50
FS-F043		Square light 220V	\$41.40		FS-F367		Respondent light 220V(by	\$87.80
FS-F044		Light 220V	\$55.30		FS-F368		Respondent light 220V(by	\$103.30
FS-F046		Round Reflector	\$8.89		FS-F040		12V 50W Mouth Lamp Bulb	\$0.65
FS-F047		Square Reflector	\$8.00		FS-F050		220V Panorama Film Viewer	\$50.30
FS-0055		Square Reflector	\$16.60		FS-F049		24V Film Viewer	\$7.87
FS-F048		220V Film Viewer	\$7.88					

Tube series




Mode	Picture	Name	FOB USD Price		Mode	Picture	Name	FOB USD Price
FS-F005		4 Hole Handpiece Tube With	\$3.81		FS-F014		4 Hole Handpiece Tube With	\$3.68
FS-F006		2 Hole Handpiece Tube With	\$3.43		FS-F015		2 Hole Handpiece Tube With	\$3.17
FS-F007		4 Hole Conneting Handpiece	\$3.81		FS-F016		Weak Suction Tube	\$1.39
FS-F008		Strong Suction Tube	\$1.78		FS-F017		3-way syringe straight tube	\$1.78
FS-F009		Tube of 2 Hole Foot Pedal	\$1.78		FS-F018		2 Hole coile Tubing(black)	\$1.65
FS-F010		Tube of 4 Hole Foot Pedal	\$3.43		FS-F019		Coile of 2 handpiece with socket	\$2.79
FS-F011		2 Hole coile Tubing	\$1.65		FS-F020		coile of triple syringe	\$1.65
FS-F012		Pipe (PU 2.5*4,3*2mm)/2	\$31.00		FS-F021		Pipe (PU 8*5mm)/10	\$50.10
FS-F013		Pipe (PU 6*4mm)/200M	\$67.80					

Valve series								
Mode	Picture	Name	FOB USD Price		Mode	Picture	Name	FOB USD Price
FS-F027		Direction Exchanger Valve	\$1.40		FS-F120		Air-Electrical Valve	\$2.66
FS-F028		Foot Valve(5m)	\$2.03		FS-F121		Water Exchange Valve(2)	\$2.54
FS-F029		Foot Valve(3m)	\$2.03		FS-F122		Strong Sucking Valve	\$2.54
FS-F030		Debris Sucking Valve	\$2.03		FS-F123		Weak Sucking Valve	\$2.54
FS-F031		Water Exchange Valve(1)	\$2.54		FS-F124		Air Regulation Valve	\$1.65
FS-F037		Simple Pressure Reducer	\$2.03		FS-F125		Water Regulation Valve	\$1.65
FS-F038		2 in 1 Air Switch	\$6.73		FS-F126		Main Air Switch	\$2.03
FS-F131		Integration Valve	\$8.77		FS-F127		Square Water Pressing	\$2.41
FS-F132		Foot Pedal Valve	\$2.41		FS-F128		Normal Open Hanging	\$1.39
FS-F133		Water Pressing Valve	\$2.92		FS-F129		Normal Close Hanging	\$1.39
FS-F134		Membrane	\$0.76		FS-F130		T-joint	\$1.39
FS-F135		Switch of Lamp	\$1.40		FS-F138		3 in 1 Valve	\$11.43
FS-F136		Air Switch	\$2.03		FS-F139		Triple Patch	\$2.41
FS-F137		Spirit Valve	\$4.19		FS-F140		4 Hole Socket	\$1.39
FS-142		4 in 1 Valve	\$22.60					

Plastic series

Mode	Picture	Name	FOB USD Price		Mode	Picture	Name	FOB USD Price
FS-F054		Holder (grey)	\$3.81		FS-F063		Rotatable Plate(round)	\$3.68
FS-F055		Holder(grey)	\$0.76		FS-F064		Rotatable plate	\$3.68
FS-F056		Holder (white)	\$1.40		FS-F065		Rotatable Plate with decoration	\$3.68
FS-F057		Air Bottle	\$1.14		FS-F066		Tissue Box	\$5.84
FS-F058		Glass	\$1.65		FS-P844		Bottle capsule	\$1.39
FS-F059		5-set Holder	\$4.95		FS-F068		Water Bottle with cover(600m	\$2.16
FS-F060		Cotton cup	\$1.52		FS-F069		Glass Frame(Luxury)	\$11.94
FS-F061		Round Protection Cover	\$3.68		FS-F070		Glass Frame	\$4.95
FS-F062		Square Protection cover	\$3.68		FS-F165		Strong Sucking Pipe	\$0.76
FS-F077		Small-small Adaptor	\$0.25		FS-F071		Weak Sucking Pipe	\$0.76
FS-F078		Big-big Adaptor	\$0.25		FS-F072		Strong Sucking Pipe	\$1.90
FS-F079		Big-big-big Adaptor	\$0.25		FS-F073		Weak Sucking Pipe	\$1.90
FS-F080		4-way Adaptor	\$0.38		FS-F074		Strong Sucking	\$0.76
FS-F081		6-way Adaptor	\$1.40		FS-F086		Weak Suction Tube(100pc	\$3.81
FS-F083		Ceramic Pipe	\$6.23		FS-F075		Arm of Lamp	\$2.66
FS-F084		Appliance Tray A	\$23.30		FS-F076		Multiple Direction	\$2.66
FS-F085		Appliance Tray B	\$23.30					













Hard series								
Mode	Picture	Name	FOB USD Price		Mode	Picture	Name	FOB USD Price
FS-LCD		LCD display	\$130.00		FS-EF03		Endoscope Frame	\$37.50
FS-EF01		Endoscope Frame	\$16.10		FS-EF04		Endoscope Frame	\$45.20
FS-EF02		Endoscope Frame	\$29.70					
FS-F051		Lamp Arm	\$38.20		FS-F141		2 Hole Socket	\$1.27
FS-F052		Column	\$13.00		FS-F142		2 Hole Cover	\$1.39
FS-F053		Balance Arm	\$62.00		FS-F143		4 Hole Cover	\$1.65
FS-F082		Set of Cuspidor	\$90.50		FS-F144		2 Hole Coupling	\$2.41
FS-F159		24V Heater	\$11.30		FS-F145		4 Hole Coupling	\$2.54
FS-F160		220V Heater	\$11.30		FS-F146		4-2 Socket	\$3.93
FS-F161		Water Filter	\$5.21		FS-F147		2-4 Socket	\$3.93
FS-F162		Dirt Collector	\$2.92		FS-F148		Straight 3 Way Syringe	\$7.75
FS-F163		Lower electric motor	\$37.50		FS-F149		3 Way Syringe	\$5.21
FS-F164		Down electric motor	\$37.50		FS-F150		Aluminium Weak Suction	\$3.68
FS-F155		Air Reducing Valve	\$6.23		FS-F151		Aluminium Strong Suction	\$3.93
FS-F156		Round Pressure Meter	\$2.03		FS-F152		Magnetic Valve	\$3.93







FS-F157		Air Reducing Valve for turbine	\$4.32		FS-F154		200W Transformer	\$15.50
FS-F158		Square pressure Meter	\$2.80					

Foot control series






Mode	Picture	Name	FOB USD Price		Mode	Picture	Name	FOB USD Price
FS-F088		Foot Control	\$28.50		FS-F181		Foot Control	\$32.27
FS-F089		Foot Control	\$28.50		FS-F182		4 Hole foot Control	\$11.68
FS-F090		Foot Control	\$23.30		FS-F183		2 Hole Foot Control	\$10.41
FS-F186		Foot Control	\$23.30		FS-F092		2 Hole Round Foot Pedal	\$8.00
FS-F091		2 hole Foot Pedal	\$5.84		FS-F093		4 Hole Round Foot Pedal	\$21.59
FS-F187		4 Hole Foot Pedal	\$11.10		FS-094		Square foot pedal	\$1.90

AIR COMPRESSOR FITTING

Mode	Picture	Name	FOB USD Price		Mode	Picture	Name	FOB USD Price
FS-F166		Pressure controller	\$6.10		FS-F173		Crash roll	\$1.39
FS-F167		Lapucitance	\$4.45		FS-F174		Round Pressure Meter(Com	\$1.90
FS-F168		220V Magnetic Valve	\$5.21		FS-F175		Rubber foot	\$0.88
FS-F169		Ball Valve	\$3.30		FS-F176		Connecting Rod(Compr essor motor)	\$11.30
FS-F170		One-way Valve	\$3.30		FS-F177		Steel Bushing(Co mpressor	\$8.25
FS-F171		Safe Valve	\$2.41		FS-F178		Bearing(Co mpressor Motor)	\$8.25

FS-AR		Air regulator	\$7.24		FS-F179		Cup leather(545 W,compress	\$4.44
FS-3W		3-way contactor	\$0.76		FS-F180		Leather(850 W,compress or Motor)	\$4.44
FS-DV		Driain valve	\$3.05		FS-IV		Intake connector	\$0.63

Other series

Mode	Picture	Name	FOB USD Price		Mode	Picture	Name	FOB USD Price
FS-F032		Light Guide(13.8mm)	\$13.30		FS-F188		Cuspidor	\$10.41
FS-F034		Curing battery	\$5.59					
FS-F035		Anti-dazzing Plate	\$0.76					
FS-F039		Bulb of Curing Machine(12V75 W)	\$2.67					